



Recombinant PD-1 [His] (DAG-WT3847)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	17.5 kDa
Format	Lyophilized
Size	1 mg
Storage	Store at -20°C

BACKGROUND

Keywords	PD-1; Programmed cell death protein 1
----------	---------------------------------------

GENE INFORMATION

UniProt ID	Q15116.3
------------	----------