



# Recombinant Human CEA protein [His] (DAG-WT3848)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Calibration, Standard
Molecular Weight	73 kDa
Format	Liquid
Size	1 mg
Buffer	20 mM PB, 300 mM NaCl, 5% Glycerol, pH 7.4.
Preservative	None
Storage	Store at -80°C

## BACKGROUND

Keywords	Human carcinoembryonic antigen; CEA
----------	-------------------------------------