



Native Human Prolactin (PRL) (>60%) (DAG-WT3849)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 60% as determined by SDS-PAGE
---------------	---------------------------------

Conjugate	N/A
------------------	-----

Applications	Calibrator
---------------------	------------

Format	Liquid or Lyophilized powder
---------------	------------------------------

Size	100 µg, 1 mg
-------------	--------------

Buffer	PBS, pH 7.4-7.6
---------------	-----------------

Preservative	0.05% sodium azide
---------------------	--------------------

Storage	Store at -20°C to -80°C
----------------	-------------------------

BACKGROUND

Keywords	Prolactin; PRL; Lactotropin; Mammotropin
-----------------	--