



Recombinant *Helicobacter pylori* Antigen [His] (DAG-WT385)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	His
Applications	ELISA, LF, TIA
Molecular Weight	60 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 µg
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C

BACKGROUND

Introduction	<p><i>Helicobacter pylori</i> is a Gram-negative, microaerophilic bacterium that can inhabit various areas of the stomach, particularly the antrum. It causes a chronic low-level inflammation of the stomach lining and is strongly linked to the development of duodenal and gastric ulcers and stomach cancer.</p>
Keywords	<p><i>Helicobacter pylori</i></p>