



Native Human Trypsin (IRT) (>25%) (DAG-WT3850)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 25% as determined by SDS-PAGE
Conjugate	N/A
Applications	Calibrator
Format	Liquid or Lyophilized powder
Size	100 µg, 1 mg
Buffer	PBS, pH 8.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Trypsin; 3.4.21.4
----------	-------------------