



Native Human HE4 (>25%) (DAG-WT3851)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 25% as determined by SDS-PAGE
Conjugate	N/A
Applications	Calibrator
Format	Liquid
Size	100 µg, 1 mg
Buffer	PBS, pH 7.4-7.6
Preservative	0.05% sodium azide
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	HE4; WFDC2; Human Epididymis Protein 4
----------	--