



Native Human Placental Lactogen (hPL) (50%) (DAG-WT3852)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	0.5
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	100 µg, 1 mg
Buffer	100 mM Tris, acetate Na, 0.5 M NaCl
Preservative	1% sodium azide
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords Human placental lactogen; hPL; Human chorionic somatomammotropin; hCS; Human chorionic somatotropin