



Native Human Serum amyloid P (SAP) (>95%) (DAG-WT3853)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity > 95% as determined by SDS-PAGE

Conjugate N/A

Applications Calibrator

Format Liquid

Size 100 µg, 1 mg

Buffer 0.3M NaCl, 20 mM Tris, pH 8.0

Preservative 0.1% sodium azide

Storage Store at 2 to 8°C

BACKGROUND

Keywords Serum amyloid P; SAP