



# Native Human Uromodulin (>90%) (DAG-WT3854)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunogen, calibrator
Format	Liquid
Size	100 µg, 1 mg
Buffer	PBS, pH 7.4-7.6
Preservative	0.05% sodium azide
Storage	Store at 2°C to 8°C for 1 year or frozen below -20°C for at least 3 years

## BACKGROUND

Keywords	Uromodulin; Tamm-Horsfall protein
----------	-----------------------------------