



## Aldosterone [BCP] (DAG-WT3861)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                |
|--------------|----------------|
| Conjugate    | BCP            |
| Applications | Immunoassays   |
| Format       | Liquid         |
| Size         | 1 mg           |
| Buffer       | PBS, pH 7.4    |
| Storage      | Store at -20°C |

### BACKGROUND

|          |                  |
|----------|------------------|
| Keywords | Aldosterone; ALD |
|----------|------------------|