



S-Adenosyl Homocysteine [BSA-Biotin] (DAG-WT3867)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA-Biotin
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4
Storage	Store at -20°C

BACKGROUND

Keywords S-Adenosyl Homocysteine; S-adenosylhomocysteine; S-adenosyl-L-homocysteine; AdoHcy