



Recombinant SARS-CoV-2 S Protein RBD (B.1.617.2) [His] (DAG-WT3869)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	34 kDa
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4
Storage	Store at -20°C

BACKGROUND

Keywords	SARS-CoV-2; S protein RBD; B.1.617.2
----------	--------------------------------------