



Recombinant *C. jejuni* Flagellin A (flaA) [His] (DAG-WT3870)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	60.4 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Campylobacter jejuni; <i>C. jejuni</i> ; Flagellin A; flaA
----------	--