



Recombinant *C. jejuni* Flagellin A (flaA) [His] (DAG-WT3873)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity > 85% as determined by SDS-PAGE

Conjugate His

Applications Immunoassays

Molecular Weight 61.5 kDa

Format Liquid or Lyophilized powder

Size 20 µg, 100 µg, 1 mg

Buffer Tris/PBS-based buffer

Storage Store at -20°C to -80°C

BACKGROUND

Keywords *Campylobacter jejuni*; *C. jejuni*; Flagellin A; flaA