



Recombinant *C. jejuni* Flagellin B (flaB) (DAG-WT3876)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity > 85% as determined by SDS-PAGE

Conjugate Unconjugated

Applications Immunoassays

Molecular Weight 59.6 kDa

Format Liquid or Lyophilized powder

Size 20 µg, 100 µg, 1 mg

Buffer Tris/PBS-based buffer

Storage Store at -20°C to -80°C

BACKGROUND

Keywords *Campylobacter jejuni*; *C. jejuni*; Flagellin B; flaB