



Recombinant *C. jejuni* subsp. *jejuni* serotype O:23/36 Uroporphyrinogen decarboxylase (hemE) (DAG-WT3881)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Purity	> 85% as determined by SDS-PAGE
Conjugate	TBD
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Campylobacter jejuni; <i>C. jejuni</i> ; Uroporphyrinogen decarboxylase; hemE
----------	---