



Recombinant *C. jejuni* subsp. *doylei* Uroporphyrinogen decarboxylase (hemE) (DAG-WT3889)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
-------------------------	---

Purity	> 85% as determined by SDS-PAGE
---------------	---------------------------------

Conjugate	TBD
------------------	-----

Applications	Immunoassays
---------------------	--------------

Format	Liquid
---------------	--------

Size	1 mg
-------------	------

Buffer	Tris/PBS-based buffer
---------------	-----------------------

Storage	Store at -20°C to -80°C
----------------	-------------------------

BACKGROUND

Keywords	Campylobacter jejuni; <i>C. jejuni</i> ; Uroporphyrinogen decarboxylase; hemE
-----------------	---