



Recombinant *C. jejuni* subsp. *jejuni* serotype O:23/36 GMP synthase [glutamine-hydrolyzing] (guaA), partial (DAG-WT3914)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Purity	> 85% as determined by SDS-PAGE
Conjugate	TBD
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Campylobacter jejuni; C. jejuni; GMP synthase [glutamine-hydrolyzing]; guaA
----------	---