



Linezolid [HSA] (DAG-WT393)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Linezolid conjugated with HSA
Species	N/A
Conjugate	HSA
Applications	Immunoassays
Format	Liquid
Concentration	1 mg/ml
Size	1 mg
Buffer	PBS
Preservative	0.05% sodium azide
Storage	Store at 2-8°C

BACKGROUND

Introduction	Linezolid is an antibiotic used for the treatment of infections caused by Gram-positive bacteria that are resistant to other antibiotics. Linezolid is active against most Gram-positive bacteria that cause disease, including streptococci, vancomycin-resistant enterococci (VRE), and methicillin-resistant Staphylococcus aureus (MRSA). The main uses are infections of the skin and pneumonia although it may be used for a variety of other infections including drug-resistant tuberculosis.
Keywords	Linezolid; LNZ; HSA