



# Recombinant *C. jejuni* subsp. *doylei* tRNA pseudouridine synthase A (truA) (DAG-WT3942)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Purity	> 85% as determined by SDS-PAGE
Conjugate	TBD
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	Campylobacter jejuni; <i>C. jejuni</i> ; tRNA pseudouridine synthase A; truA
----------	--