



Linezolid [BSA] (DAG-WT396)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|-------------------------|-------------------------------|
| Product Overview | Linezolid conjugated with BSA |
| Species | N/A |
| Conjugate | BSA |
| Applications | Immunoassays |
| Format | Liquid |
| Concentration | 1 mg/ml |
| Size | 1 mg |
| Buffer | PBS |
| Preservative | 0.05% sodium azide |
| Storage | Store at 2-8°C |

BACKGROUND

| | |
|---------------------|--|
| Introduction | Linezolid is an antibiotic used for the treatment of infections caused by Gram-positive bacteria that are resistant to other antibiotics. Linezolid is active against most Gram-positive bacteria that cause disease, including streptococci, vancomycin-resistant enterococci (VRE), and methicillin-resistant <i>Staphylococcus aureus</i> (MRSA). The main uses are infections of the skin and pneumonia although it may be used for a variety of other infections including drug-resistant tuberculosis. |
| Keywords | Linezolid; LNZ; BSA |