



# Recombinant *C. jejuni* subsp. *jejuni* serotype O:6 Peptide chain release factor 1 (prfA) (DAG-WT3987)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                  |                                                                                                                               |
|------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Product Overview | This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications. |
| Purity           | > 85% as determined by SDS-PAGE                                                                                               |
| Conjugate        | TBD                                                                                                                           |
| Applications     | Immunoassays                                                                                                                  |
| Format           | Liquid                                                                                                                        |
| Size             | 1 mg                                                                                                                          |
| Buffer           | Tris/PBS-based buffer                                                                                                         |
| Storage          | Store at -20°C to -80°C                                                                                                       |

## BACKGROUND

|          |                                                                       |
|----------|-----------------------------------------------------------------------|
| Keywords | Campylobacter jejuni; C. jejuni; Peptide chain release factor 1; prfA |
|----------|-----------------------------------------------------------------------|