



# Recombinant *C. jejuni* subsp. *jejuni* serotype O:2 Elongation factor Tu (tuf) (DAG-WT3989)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
<b>Purity</b>	> 85% as determined by SDS-PAGE
<b>Conjugate</b>	TBD
<b>Applications</b>	Immunoassays
<b>Format</b>	Liquid
<b>Size</b>	1 mg
<b>Buffer</b>	Tris/PBS-based buffer
<b>Storage</b>	Store at -20°C to -80°C

## BACKGROUND

<b>Keywords</b>	Campylobacter jejuni; <i>C. jejuni</i> ; Elongation factor Tu; tuf
-----------------	--