



Recombinant *C. jejuni* subsp. *jejuni* serotype O:2 4-hydroxy-3-methylbut-2-enyl diphosphate reductase (ispH) (DAG-WT3994)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Purity	> 85% as determined by SDS-PAGE
Conjugate	TBD
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Campylobacter jejuni</i> ; <i>C. jejuni</i> ; 4-hydroxy-3-methylbut-2-enyl diphosphate reductase; ispH
----------	---