



Recombinant Porphyromonas gingivalis Fim A [His] (DAG-WT4031)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, WB
Format	Liquid
Size	1 mg
Buffer	PBS
Storage	Store at -20°C to -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords	Porphyromonas gingivalis; Fim A
----------	---------------------------------