



Recombinant Human POLR3A protein [His] (DAG-WT4032)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	43 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris-HCl buffer
Storage	Store at -20°C to -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords DNA-directed RNA polymerase III subunit RPC1; POLR3A; RNA polymerase III subunit C1; EC:2.7.7.6

GENE INFORMATION

UniProt ID	O14802
------------	--------