



# Recombinant Human TSHR [hFc] (DAG-WT4033)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	hFc
Applications	LFIA, CLIA
Molecular Weight	109 kDa
Format	Liquid
Size	100 µg, 500 µg, 1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Keywords	TSHR; Thyroid stimulating hormone receptor
----------	--