



Recombinant Human OLR1 protein [His] (DAG-WT4038)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Reconstitution	Reconstitute at 100 µg/mL in sterile PBS
Format	Lyophilized powder
Size	50 µg
Buffer	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	LOX-1; Oxidized low-density lipoprotein receptor 1; Ox-LDL receptor 1; OLR1
----------	---

GENE INFORMATION

UniProt ID	P78380
------------	--------