



Recombinant Human SCCA1 Antigen [His] (DAG-WT4039)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	46.7 kDa
Format	Liquid
Size	100 µg, 500 µg, 1 mg
Buffer	20mM Tris, 150mM NaCl, pH8.0
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Serpin B3; SERPINB3; SCCA1; Squamous cell carcinoma antigen 1
----------	---

GENE INFORMATION

UniProt ID	P29508
------------	--------