



# Recombinant *K. pneumoniae* subsp. *pneumoniae* UDP-3-O-acyl-N-acetylglucosamine deacetylase [His] (DAG-WT4040)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                  |                                 |
|------------------|---------------------------------|
| Species          | <i>K. pneumoniae</i>            |
| Purity           | > 85% as determined by SDS-PAGE |
| Conjugate        | His                             |
| Applications     | Immunoassays                    |
| Molecular Weight | 38.1 kDa                        |
| Format           | Liquid or Lyophilized powder    |
| Size             | 20 µg, 100 µg, 1 mg             |
| Buffer           | Tris/PBS-based buffer           |
| Storage          | Store at -20°C to -80°C         |

## BACKGROUND

**Keywords** *Klebsiella pneumoniae*; *K. pneumoniae*; UDP-3-O-acyl-N-acetylglucosamine deacetylase

## GENE INFORMATION

