



# Recombinant *K. pneumoniae* OmpK36 Antigen (DAG-WT4042)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Species</b>	<i>K. pneumoniae</i>
<b>Purity</b>	> 85% as determined by SDS-PAGE
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	Immunoassays
<b>Molecular Weight</b>	38.9 kDa
<b>Format</b>	Liquid or Lyophilized powder
<b>Size</b>	20 µg, 100 µg, 1 mg
<b>Buffer</b>	Tris/PBS-based buffer
<b>Storage</b>	Store at -20°C to -80°C

## BACKGROUND

<b>Keywords</b>	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; OmpK36
-----------------	--------------------------------------------------------------

## GENE INFORMATION

<b>UniProt ID</b>	D6QLY3
-------------------	--------