



Recombinant *K. pneumoniae* Beta-lactamase SHV-3 [His] (DAG-WT4044)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	33 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Beta-lactamase SHV-3
----------	--

GENE INFORMATION

UniProt ID	P30896
------------	--------