



# Recombinant *K. pneumoniae* Undecaprenyl phosphate- $\alpha$ -4-amino-4-deoxy-L-arabinose arabinosyl transferase (a.a. 429-551) [His] (DAG-WT4045)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	17.8 kDa
Format	Liquid or Lyophilized powder
Size	20 $\mu$ g, 100 $\mu$ g, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Undecaprenyl phosphate- $\alpha$ -4-amino-4-deoxy-L-arabinose arabinosyl transferase
----------	--

## GENE INFORMATION

