



Recombinant *K. pneumoniae* mcr-8 protein (a.a. 177-565) [His] (DAG-WT4047)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	46 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; mcr-8
----------	---

GENE INFORMATION

UniProt ID	A0A2U7XX87
------------	------------