



Biotinylated Recombinant *K. pneumoniae* subsp. *pneumoniae* Chaperonin GroEL [MPB, His, Avi] (DAG-WT4048)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	MPB, His, Avi
Applications	Immunoassays
Molecular Weight	105 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Chaperonin GroEL
----------	--

GENE INFORMATION

UniProt ID	A6TH53
------------	--------