



# Recombinant *K. pneumoniae* Outer membrane porin PhoE [His, SUMO] (DAG-WT4049)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His, SUMO
Applications	Immunoassays
Molecular Weight	55 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Outer membrane porin PhoE
----------	---

## GENE INFORMATION

UniProt ID	P30704
------------	--------