



Recombinant *K. pneumoniae* Maltoporin [His, SUMO] (DAG-WT4051)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His, SUMO
Applications	Immunoassays
Molecular Weight	63.5 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords *Klebsiella pneumoniae*; *K. pneumoniae*; Maltoporin

GENE INFORMATION

UniProt ID P31242