



Recombinant *K. pneumoniae* Outer-membrane lipoprotein carrier protein [His] (Yeast) (DAG-WT4052)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	22 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Outer-membrane lipoprotein carrier protein
----------	--

GENE INFORMATION

UniProt ID	B5XYA1
------------	--------