



# Recombinant *K. pneumoniae* Outer membrane protein assembly factor BamA (a.a. 24-172) [His-B2M] (DAG-WT4053)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His-B2M
Applications	Immunoassays
Molecular Weight	33.2 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

## BACKGROUND

**Keywords** *Klebsiella pneumoniae*; *K. pneumoniae*; Outer membrane protein assembly factor BamA

## GENE INFORMATION

**UniProt ID** B5Y1J4