



Recombinant *K. pneumoniae* Outer membrane protein A (a.a. 190-344) (DAG-WT4056)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	Immunoassays
Molecular Weight	16.7 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Outer membrane protein A
----------	--

GENE INFORMATION

UniProt ID	P24017
------------	--------