



Recombinant *K. pneumoniae* subsp. *pneumoniae* Penicillin-binding protein 1B (a.a. 87-850) [His] (DAG-WT4057)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	K. pneumoniae
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	92 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Klebsiella pneumoniae; K. pneumoniae; Penicillin-binding protein 1B
-----------------	---

GENE INFORMATION

UniProt ID	A0A0H3GIH5
-------------------	------------