



Recombinant *K. pneumoniae* Fimbrial subunit type 3 [His] (Yeast) (DAG-WT4059)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	20 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Fimbrial subunit type 3
----------	---

GENE INFORMATION

UniProt ID	P12267
------------	--------