



Recombinant *K. pneumoniae* UDP-3-O-acyl-N-acetylglucosamine deacetylase [His, Myc] (DAG-WT4066)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His, Myc
Applications	Immunoassays
Molecular Weight	41.5 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; UDP-3-O-acyl-N-acetylglucosamine deacetylase
----------	--

GENE INFORMATION

UniProt ID	B5Y1U1
------------	--------