



# Recombinant *K. pneumoniae* Outer membrane protein C [His, Myc] (DAG-WT4067)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	His, Myc
Applications	Immunoassays
Molecular Weight	45 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg, 1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Outer membrane protein C
----------	--

## GENE INFORMATION

UniProt ID	Q48473
------------	--------