



# Recombinant *K. pneumoniae* Exoribonuclease 2 (partial) (DAG-WT4076)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

**Product Overview** This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.

**Species** *K. pneumoniae*

**Purity** > 85% as determined by SDS-PAGE

**Conjugate** TBD

**Applications** Immunoassays

**Format** Liquid or Lyophilized powder

**Size** 1 mg

**Buffer** Tris/PBS-based buffer

**Storage** Store at -20°C to -80°C

## BACKGROUND

**Keywords** *Klebsiella pneumoniae*; *K. pneumoniae*; Exoribonuclease 2