



# Recombinant *K. pneumoniae* subsp. *pneumoniae* Adenine phosphoribosyltransferase (apt) (DAG-WT4091)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
<b>Species</b>	<i>K. pneumoniae</i>
<b>Purity</b>	> 85% as determined by SDS-PAGE
<b>Conjugate</b>	TBD
<b>Applications</b>	Immunoassays
<b>Format</b>	Liquid or Lyophilized powder
<b>Size</b>	1 mg
<b>Buffer</b>	Tris/PBS-based buffer
<b>Storage</b>	Store at -20°C to -80°C

## BACKGROUND

**Keywords** Klebsiella pneumoniae; *K. pneumoniae*; Adenine phosphoribosyltransferase; apt