



Recombinant *K. pneumoniae* subsp. *pneumoniae* Uracil-DNA glycosylase (ung) (DAG-WT4101)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.

Species *K. pneumoniae*

Purity > 85% as determined by SDS-PAGE

Conjugate TBD

Applications Immunoassays

Format Liquid or Lyophilized powder

Size 1 mg

Buffer Tris/PBS-based buffer

Storage Store at -20°C to -80°C

BACKGROUND

Keywords *Klebsiella pneumoniae*; *K. pneumoniae*; Uracil-DNA glycosylase; ung