



Recombinant *K. pneumoniae* subsp. *pneumoniae* Succinyl-CoA synthetase subunit beta (sucC) (DAG-WT4117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Species	<i>K. pneumoniae</i>
Purity	> 85% as determined by SDS-PAGE
Conjugate	TBD
Applications	Immunoassays
Format	Liquid or Lyophilized powder
Size	1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Klebsiella pneumoniae</i> ; <i>K. pneumoniae</i> ; Succinyl-CoA synthetase subunit beta; sucC; SCS-beta
----------	--