



# Recombinant K. pneumoniae subsp. pneumoniae Glucose-6-phosphate isomerase (partial) (DAG-WT4141)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                         |                                                                                                                               |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b> | This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications. |
| <b>Species</b>          | K. pneumoniae                                                                                                                 |
| <b>Purity</b>           | > 85% as determined by SDS-PAGE                                                                                               |
| <b>Conjugate</b>        | TBD                                                                                                                           |
| <b>Applications</b>     | Immunoassays                                                                                                                  |
| <b>Format</b>           | Liquid or Lyophilized powder                                                                                                  |
| <b>Size</b>             | 1 mg                                                                                                                          |
| <b>Buffer</b>           | Tris/PBS-based buffer                                                                                                         |
| <b>Storage</b>          | Store at -20°C to -80°C                                                                                                       |

## BACKGROUND

|                 |                                                                     |
|-----------------|---------------------------------------------------------------------|
| <b>Keywords</b> | Klebsiella pneumoniae; K. pneumoniae; Glucose-6-phosphate isomerase |
|-----------------|---------------------------------------------------------------------|