



Recombinant K. pneumoniae subsp. pneumoniae Glycogen debranching enzyme (partial) (DAG-WT4151)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	This item requires custom production and lead time is between 4-12 weeks. We can custom produce based on your specifications.
Species	K. pneumoniae
Purity	> 85% as determined by SDS-PAGE
Conjugate	TBD
Applications	Immunoassays
Format	Liquid or Lyophilized powder
Size	1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	Klebsiella pneumoniae; K. pneumoniae; Glycogen debranching enzyme
-----------------	---